



STANDARD SCALP VEIN SET

STERILE, SINGLE-USE FOR INFUSION THERAPY



STANDARD SCALP VEIN SET

STERILE, SINGLE-USE FOR INFUSION THERAPY

- Ultra sharp three-beveled needle to minimize discomfort
- Siliconization of needle minimizes penetration and gliding force
- Large, soft wings for easy grip
- Kink-resistant, medical-grade PVC tubing for smooth infusion
- Color-coded to conform to ISO 6009
- Not made with latex or DEHP
- For use by healthcare professionals only



Class IIa Medical Devices

Rule 7

MDD 93/42/EEC

UMDNS 12752

PRODUCT COMPLIANCE

- CE-marked, Class IIa Medical Devices, Rule 7, MDD 93/42/ECC, UMDNS 12752
- Complies with the following norms, directives, and regulations:
 - ISO 594-1
 - ISO 594-2
 - CSN EN 1707
 - ISO 6009
 - ISO 7864
 - ISO 11607-1
 - ISO 11607-2
 - ISO 13485
 - ISO 15223-1
- Labels contain 12 languages

MANUFACTURING DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

STERILIZATION AND SHELF LIFE

Sterilization

EtO (Ethylene oxide)



Single-use only



Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that ensures sterility. This indicator comes in the form of a blue sticker that turns red when sterilized.

Shelf life

5 years



MATERIALS USED



PRODUCT RANGE OVERVIEW

| Reference code | Description | Tubing length | Color | Flow rate not less than | Priming volume |
|----------------|-------------------------------------|---------------|-------------|-------------------------|----------------|
| PSV-18-ET | Scalp Vein Set 18Gx3/4" [1.2x19 mm] | 300 mm | Pink | 46,5 ml/min | 0,46 ml |
| PSV-19-ET | Scalp Vein Set 19Gx3/4" [1.1x19 mm] | 300 mm | Orange | 36,5 ml/min | 0,47 ml |
| PSV-20-ET | Scalp Vein Set 20Gx3/4" [0.9x19 mm] | 300 mm | Yellow | 25,8 ml/min | 0,46 ml |
| PSV-21-ET | Scalp Vein Set 21Gx3/4" [0.8x19 mm] | 300 mm | Green | 20 ml/min | 0,45 ml |
| PSV-22-ET | Scalp Vein Set 22Gx3/4" [0.7x19 mm] | 300 mm | Black | 12 ml/min | 0,44 ml |
| PSV-23-ET | Scalp Vein Set 23Gx3/4" [0.6x19 mm] | 300 mm | Blue | 7,8 ml/min | 0,44 ml |
| PSV-24-ET | Scalp Vein Set 23Gx3/4" [0.6x19 mm] | 300 mm | Purple | 5,6 ml/min | 0,44 ml |
| PSV-25-ET | Scalp Vein Set 25Gx3/4" [0.5x19 mm] | 300 mm | Dark Orange | 3,7 ml/min | 0,45 ml |
| PSV-27-ET | Scalp Vein Set 27Gx3/4" [0.4x19 mm] | 300 mm | Grey | 0,35 ml/min | 0,45 ml |

PACKAGING DETAILS

| Box | Quantity (pcs) | Net weight (kg) | Gross weight (kg) | Length (mm) | Width (mm) | Height (mm) | m ³ |
|-----------|----------------|-----------------|-------------------|-------------|------------|-------------|----------------|
| Inner box | 50 | — | — | 247 | 110 | 78 | — |
| Outer box | 1000 | 5,40 | 6,40 | 572 | 263 | 344 | 0,052 |

Languages

12 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), and Finnish (FI)







Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility. Store in a dry place away from direct heat and sunlight.

Label details

| | Individual unit | Inner box | Outer box |
|---|-----------------|-----------|-----------|
| Product name | ✓ | ✓ | ✓ |
| Product size | ✓ | ✓ | ✓ |
| Quantity (pieces) | — | ✓ | ✓ |
| REF | — | ✓ | ✓ |
| LOT | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
| STERILE EO | ✓ | ✓ | ✓ |
| CE | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
| EC REP | ✓ | ✓ | ✓ |
|  | — | ✓ | ✓ |

